palajet®

Disposable Pulse Lavage System

(Disposable Electronically Pulsed Lavage Suction Apparatus)

User Manual (Instructions for Use)

Model No.: W-203

REF 5156831

STERILE EO

CE 0123

Guangzhou Clean Medical Products Manufacturing Corp.

No.163 Meidu Road, Chengjiao, Conghua District, Guangzhou, Guangdong,
China. P.C.: 510900

Importer and Distributor in EU:

Heraeus Medical GmbH

Philipp-Reis-Straße 8/13, 61273 Wehrheim, Germany

Importer and Distributor in USA:

Heraeus Medical LLC

770 Township Line Road, Yardley, PA 19067, USA

01	English	User Manual Disposable Pulse Lavage System (Disposable Electronically Pulsed Lavage Suction Apparatus) Model No.: W-203
02	French / Français	Manuel d'utilisation Disposable Electronically Pulsed Lavage Suction Apparatus N° de modèle : W-203
03	German / Deutsch	Benutzerhandbuch Elektronisch gepulstes Einweg-Lavage-Absauggerät Modell-Nr.: W-203
04	Spanish / Español	Manual del usuario Equipo desechable de aspiración y lavado mediante impulsos electrónicos N.º de modelo: W-203
05	Portuguese / Português	Manual do utilizador Disposable Electronically Pulsed Lavage Suction Apparatus Modelo n.º: W-203
06	Italian / Italiano	Manuale dell'utente Apparecchio di aspirazione per lavaggio pulsato elettronicamente monouso Nº. del modello: W-203
07	Turkish / Türkçe	Kullanıcı Kılavuzu Disposable Electronically Pulsed Lavage Suction Apparatus Model No.: W-203
08	Bulgarian / Български	Ръководство за употреба Еднократен електронно-импулсен смукателен апарат Модел №: W-203
09	Czech / Čeština	Uživatelský manuál Jednorázový systém pro pulzní laváž Modelové číslo: W-203
10	Danish / Dansk	Brugervejledning Engangsudstyr til elektronisk pulserende vaskerugning Model nr.: W-203
11	Estonian / Eesti	Kasutusjuhend Ühekordseks kasutamiseks mõeldud elektrilise impulsiga loputusvahend Mudeli nr: W-203
12	Finnish / Suomalainen	Käyttöohje Kertakäyttöinen elektronisesti pulssitoiminen huuhtelu imulaite Mallinumero: W-203
13	Greek / Ελληνικά	Εγχειρίδιο Χρήσης Συσκευή Αναρρόφησης Ηλεκτρονικής Παλμικής Πλύσης Μίας Χρήσης Αριθ. Μοντέλο: W-203

14	Hungarian / Magyar	Felhasználói kézikönyv Egyszerhasználatos elektromos jet lavage / szívókészülék Modellszám: W-203
15	Latvian / Latviešu	Lietotāja rokasgrāmata Vienreizlietojama elektroniski pulsējoša skalošanas iesūkšanas ierīce Modeļa Nr.: W-203
16	Lithuanian / Lietuvių	Vartotojo vadovas Vienkartinis Elektroninio Impulso Praplovimo Įsiurbimo Aparatas Modelio Nr.: W-203
17	Romanian / Română	Manual de utilizare Aparat de aspirare de unică folosință, cu impulsuri electronice Nr. Model: W-203
18	Slovak / Slovák	Použitý manuál Jednorazová elektrická impulzná preplachovačka Model č.: W-203
19	Slovenian / Slovenščina	Uporabniški priročnik Elektronsko impulzni sesalni aparat za enkratno uporabo Št. modela: W-203
20	Swedish / Svenska	Användarhandbok Engångsapparat med elektroniskt pulserande sug för tvätt Modellnummer: W-203
21	Dutch / Nederlands	Gebruikershandleiding Wegwerp elektronisch gepulseerd lavage-afzuigapparaat Model Nr.: W-203
22	Serbian / Srpski	Uputstvo za korišćenje Jednokratni elektronski impulsni usisni uređaj Broj modela: W -203
23	Norwegian / Norsk	Brukermanual Engangs elektronisk pulserende lavage sugapparat Modell: W-203

palajet®

CE 0123

User Manual (Instructions for Use)


STERILE EO

1. Medical Device

- 1) Device name: **palaJet® - Disposable Pulse Lavage System
(Disposable Electronically Pulsed Lavage Suction Apparatus)**
- 2) Model No.: **W-203**
- 3) **REF: 5156831**

2. Components

- 1) This device consists of **Handpiece (1)** 【Including DC motor, Compression pump, Trigger (1.1), Nozzle lock (1.2), Enclosure.】, **Nozzle (2)** 【Such as Fan Spray Nozzle (Short Nozzle) (2a), Femoral Nozzle (Long Nozzle) (2b), Femoral Brush Nozzle(Long Nozzle) (2c), Shower Spray Nozzle(Short Nozzle) (2d), Large Splash Shield(2e) or other nozzles. See the package and labels for details.】, **Waste pipe (3)** 【Including Clamp (3.1), Connecting pipe (3.2).】, **Irrigation pipe (4)** 【Including Clamp (3.1), Luer fitting (4.1).】, **Battery bag(5)**, **Power wire(6)**.
- 2) The **Vacuum extractor (7)** and **Irrigation bag (8)** are not included in this device, but prepared by medical institution. The negative pressure of **Vacuum extractor (7)** within the range of 5kPa to 40kPa shall be adjustable. The **Irrigation bag (8)** should be used with normal sterile saline or other applicable solutions (not suitable for viscous solutions).
- 3) Accessories: **Nozzles (2)**.
- 4) Detachable parts: **Nozzle (2)**.
- 5) Applied parts: **Nozzle (2)**.
- 6) Associated devices: **Vacuum extractor (7)** and **Irrigation bag (8)**.

Device structure diagram

3. Indications (Intended Use)

- 1) This device is used together with normal sterile saline or other applicable solutions, vacuum extractor. Vacuum extractor is not a must.
- 2) This device can be used for open wound cleaning, soft tissue cleaning, and surgical operation location cleaning.
- 3) This device is suitable to be used in operation room, emergency room and treatment room.

4. Working Principle

Switch the trigger to turn power on. When power is turned on, the DC motor drives the compression pump in a reciprocating motion, causing cavity volume change and discharge of the solution. The power is provided by DC batteries. There are a low power position (DC 9V), a high power position (DC 15V) and a stop tap position. Different nozzles can be fitted on according to user's need.

5. Operating Instructions

- 1) This device is portable, internal battery-powered, and consists of DC motor, compression pump, handpiece, enclosure, batteries, and nozzles.
- 2) Inspect the package. Do not use the device if package is damaged.
- 3) Consult the user manual and attached documents.
- 4) Open the package, and install **Nozzle (2)** to **Handpiece (1)**. Press down the **Nozzle lock (1.2)** to lock the **Nozzle (2)**.
- 5) Connect the **Waste pipe (3)** to the **Vacuum extractor (7)** (if prepared). Vacuum extractor is not a must.
- 6) Connect the **Luer fitting (4.1)** to the **Irrigation bag (8)**.
- 7) Hang up or place the **Battery bag (5)**. The **Battery bag (5)** should be hanged or placed on the same elevated level as the **Irrigation bag (8)** to prevent dripping of liquid from **Irrigation bag (8)** onto the battery bag.
- 8) Use the **Nozzle (2)** to aim the treatment area. Press the **Trigger (1.1)** to activate the device and start irrigation. The **Trigger (1.1)** have a low power position with symbol short “—”, a high power position with symbol long “——”, and a stop tap position with symbol “O” .

- 9) **Clamp (3.1)** can control the liquid flow in the **Irrigation pipe (4)** and **Waste pipe (3)**. Please lock the **Clamp (3.1)**, when you stop operating the device.
- 10) It can suck the waste solution concurrently while spraying or after, if **Vacuum extractor (7)** is prepared.
- 11) To press **Nozzle lock (1.2)** upward first, when withdrawing the **Nozzle (2)**.
- 12) Optional: The **Large Splash Shield (2e)** is an optional accessory to prevent splash when use **Short Nozzle**.
To install: Pinch the soft head of the short nozzle, insert into the opening of the splash shield and slightly pull in opposite directions to fixate the shield on the nozzle head.
To uninstall: Pinch the soft head of the short nozzle and remove the splash shield.

6. Specification

No.	Item	Specification
1)	Power source (Ratings)	15V DC, 10×1.5V AA (LR6) alkaline batteries.
2)	Power rating	37.5W.
3)	Handpiece dimension	About 185mm * 55mm * 40mm (Long * Wide * High).
4)	Net weight	About 0.9kg.
5)	Flow rate (high power position)	Not less than 0.9L/minute.
6)	Maximum irrigation length (high power position)	Not less than 2.5 meters.
7)	Noise (Decibel)	Less than 75dB(A).
8)	Solution form	Fog-like and thread-like.
9)	Sterilization	Sterilized using ethylene oxide before delivery.
10)	Valid period (Shelf life)	Two years.

7. IEC Classification

- 1) Internal electrical power source equipment.
- 2) Type B applied part.
- 3) IPX0.
- 4) Not category AP/APG equipment.

8. Contraindications

Not found.

9. Adverse Effect

Not found.

10. Warning and Precautions

- 1) This device should be operated only by qualified professionals.
- 2) This device should not be used after the expiry date shown on the package.
- 3) This device is for single patient use only. Do not reuse the device.
- 4) Prior to use, inspect the package. Do not use the device if package is damaged.
- 5) This device is not evaluated as category AP or APG equipment. This device is not suitable for use in the presence of flammable anesthetic mixtures with air or oxygen or nitrous oxide.
- 6) No modification of this equipment is allowed.
- 7) Except the Nozzle, there is no replaceable part and accessory within the unit.
- 8) Follow current local /national regulations to properly dispose the device after use.

11. Disposal of the Batteries

- 1) This device contains alkaline batteries. The battery will be utilizable and safe, if the device is within its useful life.
- 2) Follow local/national current regulations to properly recycle or dispose of the batteries after use.
- 3) To prevent short-circuit, Do NOT cut the power wire at any time, which may cause overheating of the battery pack and result in a potential fire risk. In case of short circuit, do not pour liquid over it nor cover it, simply leave it and dispose after it is cooled down. There is a very slim chance of short circuit if it used as instructed but proper means of disposal is highly recommended.
- 4) To safely remove the batteries, please press down on flap and remove the lid, then remove batteries from the housing.

12. Maintenance, Cleaning and Sterilization

- 1) This device is maintenance free and doesn't require routine maintenance.
- 2) This device doesn't contain user repairable parts and doesn't require repair.
- 3) This device is for single patient use only. There is no cleaning or sterilization requirement for the user. Do not resterilize.

13. Electromagnetic Compatibility Descriptions

13.1 Conform to the requirements of IEC 60601-1-2:2014 and EN 60601-1-2:2015.

13.2 Instructions for use of electromagnetic compatibility

- 1) The ME (MEDICAL ELECTRICAL) EQUIPMENT or ME SYSTEM is suitable to be used in operation room, emergency room and treatment room.
- 2) **Warning:** Don't near active HF (HIGH FREQUENCY) surgical equipment and the RF (RADIO FREQUENCY) shielded room of an ME system for magnetic resonance imaging, where the intensity of EM (ELECTROMAGNETIC) disturbances is high.
- 3) **Warning:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- 4) This device doesn't contain user repairable cables, transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION.
- 5) **Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Disposable Electronically Pulsed Lavage Suction Apparatus

(W-203), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

- 6) **NOTE:** The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

13.3 Technical description of electromagnetic compatibility

- 1) All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to ELECTROMAGNETIC DISTURBANCES for the EXPECTED SERVICE LIFE. Disposable use product. Do not change any parts of the product.
- 2) In the table below the guidance and manufacturer's declaration -electromagnetic emissions and immunity are shown:

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions (CISPR 11)	Group 1
RF emissions (CISPR 11)	Class A
Harmonic emissions (IEC 61000-3-2)	N/A
Voltage fluctuations/ flicker emissions (IEC 61000-3-3)	N/A

Table 2

Guidance and manufacturer's declaration - electromagnetic immunity		
Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 6 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV
Electrical fast transient / burst IEC 61000-4-4	Power supply lines: ± 2 kV input/ output lines: ± 1 kV 100 kHz repetition frequency	N/A
Surge IEC 61000-4-5	line(s) to line(s): ± 0.5 kV, ± 1 kV. line(s) to ground: ± 0.5 kV, ± 1 kV, ± 2 kV.	N/A
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_T 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T 1 cycle And 70% U_T 25/30 cycles Single phase: at 0 0% U_T 250 cycle	N/A
Power frequency magnetic field IEC 61000-4-8	50Hz: 30 A/m 60Hz: 30 A/m	50Hz: 30 A/m 60Hz: 30 A/m
Conducted RF IEC 61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	N/A
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2700 MHz 80 % AM at 1 kHz	3 V/m 80 MHz to 2700 MHz 80 % AM at 1 kHz

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Table 3

Guidance and manufacturer's declaration - electromagnetic immunity							
	Test frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	Immunity test level (V/m)
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
	450	430– 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
1970							

	2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
	5500						
	5785						

14. Environmental Requirements

1) Storage conditions:

This device should be stored in non-corrosive gases, cool, dry, and ventilated place.

Temperature: -10°C~+40°C; Humidity: ≤95%RH; Air pressure: 50kPa~106kPa.

2) Operating conditions:

Temperature: 5°C~30°C; Humidity: ≤80%RH; Air pressure: 86kPa~106kPa.

15. Symbol Descriptions

The following symbols will appear on user manual, labels and packages. Some of the symbols represent standards and compliances associated with the device and its use.

Sterilized using ethylene oxide	Do not re-use	Consult instructions for use
Do not re-sterilize	Caution	Do not use if package is damaged and consult instructions for use
CE marking of conformity	Authorized Representative in the European Community/ European Union	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
The law restricts this device to sale, distribution, and use by or on the order of a physician.		
Model number	Catalogue number	Unique device identifier
Batch code (Lot number)	Date of manufacture	Use-by date
Manufacturer	Medical device	Type B applied part
This symbol on the trigger side means stop tap position	Importer	Direct current
This symbol on the trigger side means low power position (DC 9V)	Distributor	Temperature limit
This symbol on the trigger side means high power position (DC 15V)	This way up	Fragile, handle with care
UK Responsible Person	Keep dry	Keep away from sunlight

16. Manufacturer

Name: **Guangzhou Clean Medical Products Manufacturing Corp.**

Address: No.163 Meidu Road, Chengjiao, Conghua District, Guangzhou, Guangdong, China.

P.C.: 510900

Tel: +86-20-87654060, 87817888

Fax: +86-20-87303882

E-mail: info@cleanmfg.com

Web: www.cleanmfg.com

17. EU Authorized Representative

Name: **SUNGO Europe B.V.**

Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands.

Tel: +31 (0)20 2111 106

E-mail: ec.rep@sungogroup.com

18. Importer and Distributor in EU

Name: **Heraeus Medical GmbH**

Address: Philipp-Reis-Straße 8/13, 61273 Wehrheim, Germany.

19. Importer and Distributor in USA

Name: **Heraeus Medical LLC**

Address: 770 Township Line Road, Yardley, PA 19067, USA.

20. UK Responsible Person

Name: **SUNGO Certification Company Limited**

Address: 3rd Floor, 70 Gracechurch Street, London, EC3V 0HR, England, United Kingdom.

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Navn: **Guangzhou Clean Medical Products Manufacturing Corp.**
Navn: **SUNGO Europe B.V.**
Navn: **Heraeus Medical GmbH**